

Low Dose Radiographic System for Pediatrics

Y. KOVALENKO¹, S. MIROSHNICHENKO³, A.NEVGASIMY⁴, S.BALASHOV³, N.POTRAHOV²,
Y.MIRONOVA⁵

¹National Medical Academy of Postgraduate Education named after P.L. Shupyk, Kiev, Ukraine

²St. Petersburg Electro-Technical University, St. Petersburg, Russia

³National Aviation University, Kyiv, Ukraine

⁴Teleoptic PRA, Kiev, Ukraine

⁵KRU KTMO "University Hospital", Simferopol, Ukraine

Abstract — The paper deals with a specialized digital X-ray system for use in children's trauma centers and maternity hospitals. The basic requirements for such equipment: high contrast and spatial resolution of the images with minimal radiation exposure to the baby. It is proposed the set of equipment which meets these requirements: the microfocus X-ray generator (output power of less than 50 W, the focal spot of less than 50 microns), a special tripod and a digital receiver on a array of the photodiode matrices, which allows to obtain a digital image with a spatial resolution to 9.0 line pairs per millimeter with entrance exposure dose in the plane of the receiver is less than 500 mR

Index Terms — pediatrics, low-dose radiographic system, high contrast and spatial resolution, microdose radiation exposure.

I. PURPOSE

In pediatric radiology patients are significantly different from each other both by means of size, and on the density structure. On the one hand, young children have not formed bone structure, geometry objects of study are small, which requires a high contrast and spatial resolution of the images. On the other hand, the younger the child, the more sensitive it's organism to radiation. This makes ultrasonic method a preferred application when examining young children. However, in many cases, X-ray method should be used, & it is necessary to minimize the effects of radiation to the patient. The purpose of the work was to create a specialized X-ray equipment to obtain diagnostic images with high contrast and spatial resolution at microdose (several hundreds of micro-roentgen) radiation exposure to the patient.

II/ MATERIALS AND METHODS

Digital X-ray system consists of a small mobile universal tripod with a special table-stand, a micro transmitter and a highly sensitive digital receiver based on photodiode matrices array.

By changing the position of a radiolucent tabletop between the transmitter and the receiver X-ray studies of objects in X-ray mode, zoom in both the horizontal and vertical planes could be performed (Fig.1a,b).




a

b


Fig. 1. Low-dose X-ray system:
a – in horizontal position, b- in vertical position

III. RESULTS


Small-sized mobile radiographic system that provides high-quality digital images of X-ray generator with a power less than 50 Wt was implemented. The subject can be located directly at the receiver and away from it.



a



b



c

Fig. 2 Panoramic X-ray image of mammographic phantom (a) and X-ray images of some of its parts made in magnification mode (b,c)

By changing the distance from the radiation source to the examined object (f_1) and the distance from the object to the receiver (f_2), you can change the coefficient of magnification $C_m = (f_1 + f_2) / f_1$. Use the zoom mode allows you to change the spatial resolution of the images. The larger the zoom, the higher the spatial resolution, but smaller field of view. For illustration, Figure 2 shows a snapshot of mammographic phantom, and the 2b - detailed pictures of the individual parts, made in zoom mode.

At the X-ray images of mira images with a spatial resolution of up to 9.0 line pairs per millimeter were obtained. In this case, the exposure dose in the plane of the receiver does not exceed 500 mR (microdose mode). To increase the sensitivity (lower dose) there is a binning mode in the digital receiver, where the signal of several pixels are merged into one. Binning modes and zoom extend the operator's opportunities to choose the optimal ratio "image quality-dose", which is especially important in significant variability of research objects. The possibility of a detailed analysis of the integrity of bone tissue using microfocus X-ray diffraction in combination with the assessment of soft tissue injuries using ultrasound allows to make more accurate diagnosis of traumatic injuries. A low radiation load on the object of study are responsible for the feasibility of using Low dose x-ray systems to take the X-ray images of babies who have increased sensitivity to ionizing radiation.

IV. CONCLUSION

Application of considered digital microfocus radiographic system is most promising in hospitals for x-ray examinations in infants and in children's emergency station in combination with ultrasound scanners.

Digital Radiographic Complexes for Traumatological Department

V.SUDAKEVICH¹, Y.KOVALENKO², S.MIROSHNICHENKO³, S. BALASHOV³, A. NEVGASIMY⁴

¹Uman central regional hospital

²Natsionalnaya Medical Academy of Postgraduate Education named after P.L. Shupyk, Kiev, Ukraine

³National Aviation University, Kyiv, Ukraine

⁴Teleoptic PRA, Kiev, Ukraine icnbme@mail.utm.md

Abstract — The paper deals with a budget digital X-ray equipment in operating trauma unit. The proposed set of equipment lets receive X-ray images of required quality at a significantly reduced requirements for radiation protection in the operating room, and to protect the personnel. We consider two variants of the equipment: for operating rooms in the traumatological department in the central district hospitals on the basis of portable X-ray unit and pediatric trauma unit on the basis of micro focus X-ray generator.

Index Terms — *casualty departments, high-frequency X-ray monoblock, microfocus X-ray unit, digital detector, minimal price, radiation safety*

I. PURPOSE

Radiological imaging plays a critical role in patients with trauma. On the one hand, it allows to diagnose traumatic injuries, on the other hand, to control the quality of surgical operations. Currently, expensive X-ray equipment with C-arc is widely used, but it is a significant radiation exposure to patients and staff. The aim of the work is to develop and implement to the clinical practice the budget options of digital X-ray equipment for rapid diagnostics and monitoring performance in quasi-real time the results of surgical interventions in casualty departments with minimal radiation exposure to patients and staff.

II. MATERIALS AND METHODS

The kits consisting of: high-frequency X-ray monoblock SY-HF-110 (output power 2.2 kW) with a mobile stand, digital detector "Iona-R4000" with stand S-30TS and adapter to the surgical table for operating (Fig.1) and microfocus X-ray unit type "Pardus", a special tripod and digital detector "Iona-R2430" for pediatric & wrist trauma departments were offered (Fig. 2).



Fig.1 Digital X-ray complex in the operating room of casualty department



Fig.2 Microfocus digital X-ray complex

Both complexes can carry out examinations in direct and lateral views. Microfocus X-ray generator allows to perform X-ray images in zoom mode, and thereby enhance its detail.